


WIRKPRINZIP ÜBER DEN MUND METHOD OF ACTION THROUGH THE MOUTH

ANWENDUNG / USE

RC-FIT® DUO ÜBER DEN MUND

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RC-FIT® DUO THROUGH THE MOUTH

A) TRAINING & THERAPIE WITH VIBRATIONS

As it is in sports and everyday life, having more oxygens means having a greater ability to function.
The Inlay with integrated valve produces oscillations that are acoustically perceived as sounds. These spread in the airways. As a result, shortness of breath is decreased during physical exertion. Mucilage and coughing in acute and chronically ill patients can be reduced. Through stabilisation of the airways more oxygen can be taken in. In addition, body function perception is achieved.
The goal is to support the cleaning function of the airways in the best possible manner. This can be very effective within the scope of simple colds all the way to chronic pulmonary diseases. Aside from moistening the airways with fresh, moist air (if there is no fever present), inhalations, or moistening nasal sprays in order to keep the cleaning mechanism working, the RC-FIT® therapy device can also free the airways of stubborn mucus. The breeding ground for infections is thus minimised.

B) ANCHOR FOR SHORTNESS OF BREATH

As with the pursed lips technique, the RC-FIT® DUO through the mouth can stabilise the bronchi in cases of shortness of breath, through a deliberately slow exhalation. Hyperventilation can be reduced or prevented. Thanks to its practical pocket shape, the device is always close at hand in the lowest intensity level.
Practice removing with the RC-FIT® DUO through the mouth and the use of pursed lip technique about 60 times to be optimal for emergencies to be armed.

USE

Use of the RC-FIT® DUO through the mouth is not dependent on the force of gravity. The device may be used in a reclining, sitting or standing position. For the best effect, however, using the device while standing is preferable. Drink a glass of water before training in order to keep your airways moist.
I. Please inspect the device for foreign bodies or damages prior to every use. If any parts are missing or damaged, replace them immediately.
II. Always start your first training session in the lowest intensity level "level I" (largest-possible triangle).
The smaller the triangle becomes, the higher your training and therapy load (intensity). RC-FIT® DUO through the mouth has a total of 4 intensity levels.

Average pressure during exhalation / inhalation

Intensity level	Flow rate (pressure difference to ambient pressure)
I	5-9 l/min
II	10-19 l/min
III	20-29 l/min
IV	30-39 l/min

Frequency during inspiration: 165 Hz +4%
Frequency during exhalation: 130 Hz +4%

III. A) Place the RC-FIT® DUO through the mouth in your mouth and close your lips around the mouthpiece. Please do not close the triangle with your fingers.

Breathe in and out through the mouth.

If the sound at the start of each breathing manoeuvre is loud and then quickly becomes quieter, please adjust your breathing accordingly.

Make sure that the sound when breathing through the RC-FIT® DUO through the mouth is uniform and continuous. In such a case you will know that you have surmounted the training load.
Concentrate on the sounds occurring with your breathing. Try to produce these sounds to be as balanced as possible, for as long as possible.

Depending on your training condition, take a 1-minute break after 1½ minutes. One training unit should consist of 10–30 breaths, 3 x daily. In the case of Parkinson's with dysphagia, training of 5 x daily with 5 exhalations each time should be conducted.

B) To train for cases of shortness of breath or while experiencing shortness of breath always use only the "I" intensity level.
Try to breathe in the goalkeeper position or the cat diver position calmly and in a prolonged manner, via the RC-FIT® DUO through the mouth. Repeat this breathing cycle approx. 10 times if possible.

If the training is too light, or the sound is too loud even after 30 breaths, increase the intensity (resistance) by rotating the mouthpiece, see Fig. II.

CLEANING AND DISINFECTION (Boiling)

If used regularly, clean and disinfect the device once a week, please clean daily if there is a lot of mucus. If necessary (e.g. if no more sounds are produced), insert parts of the RC-FIT® DUO through the mouth in a warm, light detergent solution and then boil out.

1. Cleaning and disinfection
For cleaning purposes, pull off the mouthpiece from the intensity regulator. Then, carefully remove the inlay with integrated valves from the intensity regulator.

IMPORTANT! Try to avoid any contact with the integrated valves! Wash all parts in a solution of lukewarm clean water and standard liquid soap detergent. Rinse with clear water.

After this cleaning, disinfect all parts in a vaporisator or in boiling water for approx. 5 minutes.

2. Drying
In order to prevent mould formation, all parts must be completely dried, to accelerate this process, shake off any excess water before setting the parts down. Allow all parts to air dry. Do not rub dry the inlay with towels.

3. Assembly
Before assembly, please wash your hands and ensure that all parts of the RC-FIT® DUO are clean. There may be no dust residue or scattered particle deposits, as these could be breathed in and reach the lungs.

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Recommended therapy intervals for the respiration:
RC-FIT® DUO 3 x daily 10–30 breaths, 7 days a week. Additional usage when you notice heavy mucus obstruction or occasional shortness of breath.

TIPPS:
► Integrate the therapy and training into your daily routine, especially when you notice very little mucus, e.g. always before brushing your teeth.
► To obtain a better effect from your medication spray, integrate this inhalation into therapy with the RC-FIT® DUO.

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RC-FIT® DUO THROUGH THE NOSE

POSSIBLE EFFECTS

- Reactivates natural breathing
- Improves ventilation of the sinuses and middle ear
- Reduces shortness of breath, caused by e.g. psychological factors
- Reduces the breathing rate
- Strengthens the respiratory muscles
- Improves gas exchange
- Frees the upper airways from irritating phlegm
- Reduces coughing due to post-nasal-drip
- Reduces bronchial hyperactivity at asthma
- Enhances drug deposition

USE

Use of the RC-FIT® DUO through the nose is not dependent on the force of gravity. The device may thus be used in a reclining, sitting and standing position. According to the guidelines*, we recommend cleaning your nose with a nasal rinse before using the RC-FIT® DUO through the nose if you have sinusitis. Alternatively, you can use a moisturising nasal spray (e.g. RC-Pflege® N). Experience has shown that this results in a faster and more pronounced therapeutic effect.
I. Please inspect the device for foreign bodies or damages prior to every use. If any parts are missing or damaged, replace them immediately.
II. Always start your first training session in the lowest intensity level "MINIMUM" (largest-possible triangle).

The smaller the triangle, the higher your training and therapy load (intensity). The intensity levels on the nosepiece are continuously adjustable.

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III. Please clean your nose before beginning the training. Hold the RC-FIT® DUO through the nose as shown. Insert the round opening of the nosepiece of the RC-FIT® DUO through the nose into one nostril and close the other nostril with the index finger of your other hand. Ensure that the triangle of the intensity control remains open at all times.

Breathe in through your nose and after a short break, breathe out through your nose again.

Repeat the process with the other nostril after completing the training session (approx. 10–30 breaths).

If the sound at the start of each breathing manoeuvre is loud and then quickly becomes quieter, please adjust your breathing accordingly. Make sure that the sound is uniform and continuous when breathing through the RC-FIT® DUO through the nose. By this, you will know that you have surmounted the training load. Concentrate on the sounds created by your breathing. Try to produce these sounds as balanced as possible, for as long as possible. The sounds may be slightly delayed depending on the intensity. If no sound is heard, try adjusting the angle slightly and start the breathing very softly.

Depending on your training condition, take a 1-minute break after 1½ minutes. One training unit should consist of 10–30 breaths, 3 x daily, per nostril.

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If used regularly, clean and disinfect the device once a week, please clean daily if there is a lot of mucus. If necessary (e.g. if no more sounds are produced), insert parts of the RC-FIT® DUO through the nose in a warm, light detergent solution and then boil out.

1. Cleaning and disinfection
For cleaning purposes, pull off the nosepiece from the intensity regulator. Then, carefully remove the inlay with integrated valve from the intensity regulator.

IMPORTANT! Try to avoid any contact with the integrated valve! Wash all parts in a solution of lukewarm clean water and standard liquid soap detergent. Rinse with clear water.

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Recommended therapy intervals for the respiration:
RC-FIT® DUO 3 x daily 10–30 breaths, 7 days a week. Additional applications whenever needed.

TIPPS:
► Integrate the therapy and training into your daily routine, especially when you notice very little mucus, e.g. always before brushing your teeth.

* Leitlinie der German Otorhinolaryngology Society, Head and Neck Surgery and the German Society of General and Family Medicine, DGAH (online) (Naseninnenhöhlenentzündung)