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1. Review by Veasey SC, Guilleminault C, Strohl KP et al. Medical Therapy for Obstructive Sleep Apnea: A Review by the Medical Therapy for Obstructive Sleep Apnea Task Force of the Standards of Practice Committee of the American Academy of Sleep Medicine. *Sleep* 2006; 29(8):1036-1044.
2. Penzel T, Möller M, Becker HF, et al. Effect of sleep position and sleep stage on the collapsibility of the upper airways in patients with sleep apnea. *Sleep* 2001; 24:90-95.
3. Oksenberg A, Khamaysi Y, Silverberg DS, et al. Association of body position with severity of apneic events in patients with severe nonpositional obstructive sleep apnea. *Chest* 2000; 118:1018-1024.
4. Sahlin C, Svanborg E, Stenlund H, et al. Cheyne-Stokes respiration and supine dependency. *Eur Respir J* 2005; 25: 829-833.
5. Cartwright RD. Effect of sleep position on sleep apnea severity. *Sleep* 1984; 7: 110-14.
6. Oksenberg A, Silverberg DS, Arons E, et al. Positional vs nonpositional obstructive sleep apnea patients. Anthropomorphic, nocturnal polysomnographic and multiple sleep latency test data. *Chest* 1997; 112:629-639.
7. Mador MJ, Kufel TJ, Magalang UJ, et al. Prevalence of positional sleep apnea in patients undergoing polysomnography. *Chest* 2005; 128:2130-2137.
8. Nakano H, Ikeda T, Hayashi M, et al. Effects of Body Position on Snoring in Apneic and Nonapneic Snorers. *Sleep* 2003; 2:169-172.
9. Skinner MA, Kingshott RN, Filsell S, et al. Efficacy of the 'tennis ball technique' versus nCPAP in the management of position-dependent obstructive sleep apnoea syndrome. *Respirology* 2008; 13: 708-715.
10. Epstein LJ, Kristo D, Strollo PJ, et al. Clinical guideline for the evaluation, management and long-term care of obstructive sleep apnea in adults. *J Clin Sleep Med* 2009; 5(3):263-276.
11. Hidalgo Armas L, Turino C, Cordero-Guevara J et al. A new postural device for the treatment of positional obstructive sleep apnea. A pilot study. *Respir Med* 2019; 151:111-117.

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T. 933 90 50 30

www.somnibel.eu · somnibel@sibelmed.com

This manual can be purchased through the After Sales Service.

COPYRIGHT

No part of this publication may be reproduced, transmitted, transcribed, stored in a back-up system or translated into any language or computer language in any form or by any means, electronic, mechanical, optical, chemical or manual without the express written consent from SIBEL S.A.U.

DISCLAIMER

SIBEL S.A.U. is responsible for the security, reliability and performance of this equipment only if:

- The place where the system is installed or used meets the requirements for electrical installations IEC and other applicable regulations.

C €0197

PRODUCT IN COMPLIANCE WITH MEDICAL DEVICE DIRECTIVE 93/42/EEC
(CLASS IIa).

Thank you for choosing **Sornibel**!

This product has been designed by the RDI department of SIBEL S.A.U. in cooperation with the sleep units of the University Hospital Araba / Osakidetza (Vitoria, SPAIN) and of the University Hospital Arnau de Vilanova (Lleida, SPAIN). **Sornibel** is manufactured according to the European Medical Device Directive and with the best quality standards.

Revised
Date: 2019-09
Technical Director

Approved
Date: 2019-09
Sales Director

INDEX

1. SAFETY	4
1.1 INTENDED USE	4
1.2 INDICATIONS FOR USE	4
1.3  WARNINGS.....	6
1.4 DISPOSAL OF ELECTRICAL AND ELECTRONIC DEVICES BY DOMESTIC USERS IN THE EUROPEAN UNION.....	9
2. INSTALLATION INSTRUCTIONS	10
2.1 MAIN FEATURES	10
2.2 PACKING LIST	10
3. INSTRUCTIONS FOR USE	11
3.1 SOMNIBEL FUNCTIONING.....	12
3.2 COMPLIANCE ENHANCEMENT FUNCTION	13
3.3 VISUAL INDICATIONS AND WARNINGS.....	14
3.4 TRAVELLING WITH THE SOMNIBEL	15
4. CLEANING INSTRUCTIONS	15
5. TECHNICAL SPECIFICATIONS	16
6. APPLICABLE DIRECTIVES AND STANDARDS	18
7. SYMBOLS.....	19
8. BIBLIOGRAPHY	20
Annex 1: ELECTROMAGNETIC COMPATIBILITY	21

1. SAFETY

Somnibel is a new system for the treatment of patients with Sleep-Related Breathing Disorders (SRBD), mainly the Sleep Apnea-Hypopnea Syndrome (SAHS) or habitual snoring during sleep. It has been developed by the RDI Department of SIBEL S.A.U. with the collaboration of reference Health Centers and doctors specialized in Sleep Disorders.

Somnibel has been designed for maximum safety. The complete instructions for use should be read before proceeding with the device use. Failure to do so may result in injuries to the user and damage to the device and/or accessories.

1.1 INTENDED USE

The **Somnibel** device is intended for the treatment of positional Sleep Apnea Hypopnea Syndrome (SAHS) or positional snoring by applying a gentle vibrating stimulus when the patient sleeps in the supine position to induce a change in the body position, thus reducing the incidence of sleep respiratory events.

The following conditions must be taken into account:

- Body-worn device
- Use into the patient's home or similar indoor use (not for outdoor use).
- Not intended for use in moving transport vehicles

1.2 INDICATIONS FOR USE

It is well known that more than 50% of SAHS patients experience a worsening in the frequency and severity of respiratory events when sleeping in the supine position, compared to other body positions^{1, 2, 3, 4}. This positional dependency is more prevalent in mild to moderate, younger and less obese patients, which are the majority of the total population of SAHS patients^{5, 6, 7}. Positional SAHS is defined as an abnormal Apnea-Hypnea Index (AHI) that is double in supine position compared to other positions^{5, 7}. Additionally, the severity and frequency of snoring is also positional dependant in the majority of habitual snoring subjects⁸. When the supine position is avoided, the frequency of respiratory events is reduced in patients with positional Obstructive Sleep Apnea (OSA)⁹ and positional therapy is considered an effective treatment for OSA in selected patients¹⁰.

Somnibel is a body-worn, active device able to detect the position of the patient during sleep and to generate a gentle vibration, depending on the detected position, that induces the patient to change his/her position.

Somnibel has been designed to be used by the patient at home while lying down to sleep and can be placed at the forehead or chest. In both placements, the device vibration is able to induce a position change and, consequently, a reduction of respiratory events. **Somnibel** includes several compliance enhancement features for favoring the patient's sleep and making it more comfortable and easy to use. The use of Somnibel has proven to be effective in reducing the AHI and the time in supine position and in improving the oxygen saturation, without affecting the quantity nor quality of sleep in patients with positional OSA¹¹.

Somnibel is indicated for patients older than 18 years, weighing over 40 kg and a height over of 145 cm. Although not required, it is recommended to use the device under the supervision of a sleep physician for obtaining better treatment results.

Somnibel uses a rechargeable Lithium-polymer battery and includes USB connection for battery charge. It is recommended to fully charge the battery before the first use of the device.

Somnibel must be adhered to the patient's body, at the forehead or chest, by means of a specifically designed disposable fastening adhesive. The fastening adhesive should not be placed onto injured skin. The patient should adopt a lying sleeping position (trunk inclination $< 30^\circ$) while using the device.

Although **Somnibel** can be used within conventional environmental conditions (see section TECHNICAL SPECIFICATIONS), it is recommended to sleep with a room temperature between 18°C and 26°C . A hot environment may induce excessive sweating that can deteriorate the adherence of the adhesive.

The use of the device in an environment with vibrations, such as moving transport vehicles (car, train, airplane, etc.) may modify its functioning.

1.3 WARNINGS

Somnibel is not intended to be used with other energy sources that are not covered in this manual.

Use only **Somnibel** with accessories provided by the manufacturer or dealer, or those that meet the specifications of this manual. The use of other accessories with **Somnibel** can cause damage to the device or affect the measurements.

The battery charger connector includes a protecting cover for avoiding penetration of liquids into the device and for avoiding electrostatic discharges during their use. Do not use this accessory without this cover.

Contact of liquids with the internal parts of the device and the connectors must always be avoided. The system is resistant to moderate splashing and dripping (Protection level IP22: protected against solid objects with a diameter of 12.5 mm and above; protected against water drops falling vertically with a maximum inclination of 15 degrees).

Accessories should be handled by their strongest parts, which are the connectors. They should also not get wet or exposed to very abrupt changes of temperature. Do not apply excessive stress to the accessories. In particular, avoid pulling or bending any part of the cables.

Do not disconnect the accessories from the device by pulling the cable. You can damage the product or accessories reducing the product's safety. Hold always the connectors when disconnecting the device.

Deteriorated accessories (battery charger with broken case, cable or connector) should not be used since there is risk of cramps. Contact SIBEL S.A.U. or your provider for acquiring new ones.

Deteriorated fastening adhesives (excessively wrinkled or partially removed from the protecting cover) should not be used since their adhesiveness may be reduced. Dispose defective adhesives and use a new one.

In case of device overheating during battery charge or during use, immediately disconnect it from the mains and stop using the device. Contact SIBEL's After Sales Service or your dealer.

Fastening adhesives are intended for single-use. Do not reuse the adhesives since their adhesiveness is considerably reduced after the first use. Use always a new adhesive for each treatment session.

The cleaning instructions in this manual must be carefully followed.

No parts are allowed for temporary immersion. IT MAY CAUSE ELECTRIC DISCHARGE.

The detection of body position and the generation of vibrating stimulus are considered essential operation (so that the specified accuracy can be met), as well as the attachment of the device to the patient's body and the visual indications.

The **Somnibel** device is prepared to work at room temperature. Avoid exposing any part of the device or accessories to heat sources. Also avoid direct exposure to sunlight. Temperature changes cause condensation and moisture. Before using the device, allow the it to acclimate to ambient temperature.

The device is not designed to work in an explosive environment or in the presence of flammable anesthetics or gases of any kind. **MAY CAUSE EXPLOSION.**

The device should be used in an acoustic environment that allows the patient to sleep normally.

Keep your device protected from shock and vibration. During transportation, place all the items in the original packaging. The material provides enough protection against small accidental impact.

Do not throw the packaging after using the device. It can be used for device transportation.

Do not use **Somnibel** or the fastening adhesive continuously for a period longer than 10h. Prolonged use of the product can induce adverse effects such as headache, skin marks or redness after adhesive removal. If you detect any of these effects, stop using the device and contact your physician.

Do not try to open the device. In case of malfunction or unexpected operation, please contact SIBEL's After Sales Service or your dealer.

In case of receiving **Somnibel** with a deteriorated packaging, do not use the device and contact your courier agency or SIBEL's After Sales Service.

Do not use the device in and MRI environment.

The use of mobile phones, transmitters and similar equipment generating radio frequency emissions and placed next to the device is not allowed during treatment with **Somnibel**. Therefore, do not use the device in the presence of radio equipment (e.g. mobile phones). Follow the recommendations regarding the separation distance specified in the manufacturer's declaration for EMC in this manual (Annex 1).

The use of **Somnibel** in combination with a pacemaker is not recommended. Consult your medical specialist before using **Somnibel**.

1.4 DISPOSAL OF ELECTRICAL AND ELECTRONIC DEVICES BY DOMESTIC USERS IN THE EUROPEAN UNION

Never dispose **Somnibel** in the household trash, since it contains a Lithium-Polymer battery and electronic components. It must be disposed properly and may need to be recycled in accordance with the statutory requirements in your country.

- Waste of batteries Directive: the lithium-polymer battery contained into the device must be extracted and disposed appropriately before disposing the device.
- Materials according to the RoHS Directive: the device and its accessories are RoHS compliant.
- Materials according to the Medical Device Directive: none of the components or accessories that may be used in combination with the device contain phthalates. The device and all the accessories are latex free.
- Materials according to the REACH regulation: neither the device nor its accessories use any hazardous substance according to REACH regulation.

Fastening adhesives for **Somnibel** can be disposed in the household trash, unless local legal requirements are applicable. See the fastening adhesive User's Manual for additional information on disposal. Information on proper disposal Information on proper disposal is available from SIBEL's After Sales Service or your dealer.

2. INSTALLATION INSTRUCTIONS

2.1 MAIN FEATURES

- Automatic ON/OFF
- Internal rechargeable battery: > 3 nights (10h/night)
- Visual indicator of battery status
- Initial self-test

2.2 PACKING LIST

CODE 08767 QTY. 1 SOMNIBEL POSITIONAL THERAPY DEVICE

CODE 08741 o 08742 QTY. 1 USB CHARGER

CODE 08741 TRUMPower TMW7-5-E-UB
CODE 08742 GLOBTEK GTM41134-0606-1.0

CODE 08753 QTY. 1 30u. FASTENING ADHESIVES BOX, includes:

- 30 FASTENING ADHESIVE

- ADHESIVE USER'S MANUAL



CODE 08754 QTY. 1 Sominibel USER'S MANUAL

3. INSTRUCTIONS FOR USE

1 Clip one fastening adhesive to a **Somnibel** device until you hear a "click" (the light blinks 3 times in green and the device vibrates).

3 Remove the protective cover from the adhesive.

2 Pay attention to the indications of the device after turn ON (consult section VISUAL INDICATIONS AND WARNINGS).

4 Pay attention to the correct orientation indicated by an arrow at the label of the device. Place the device and adhesive to the forehead or chest.

5 Lie down for sleeping.

6 On the next morning, remove the fastening adhesive from the skin and unclip it from the device (check that the green light blinks twice).

7 Dispose the fastening adhesive and its cover.

8 Connect the device to the mains by means of the battery charger. Check that the light turns on green.

9 When the light turns off, the device battery is fully charged and it can be disconnected from mains.

10 When not in use, store the device and accessories into its packaging.

Do not forget to unclip the fastening adhesive from **Somniael** after each use. This will expand battery life.

If **Somniael** indicates low battery after clipping the fastening adhesive, connect the device to mains before using it for assuring a complete therapy during the whole night.

3.1 SOMNIBEL FUNCTIONING

Somniael is automatically turned ON when the fastening adhesive is clipped to the device. At the beginning, the device performs a self-test for checking its correct functioning. The visual indicator on the device will provide indications of the results of self-test and battery status (see sections VISUAL INDICATIONS AND WARNINGS).

After turning ON, **Somniael** will wait for 15 minutes before starting the stimulus application so that the patient can restfully get asleep.

Somniael device generates a gentle vibrating stimulus when the supine position is detected for inducing the patient to change his/her body position. If the patient does not change the position, the device will wait for 30 seconds and then it will increase the stimulus intensity. The intensity will be progressively increased until the patient changes the position (4 intensity levels).


3.2 COMPLIANCE ENHANCEMENT SETTINGS

The **Somnibel** device includes additional features for enhancing patient compliance to the positional therapy.

These features are configured with default parameters, as follows:

Settings	Default value	Description	Enhancement
EasySleep Delay	15 min	Time between device ON and activation of stimulus if patient is in supine position.	Allows the patient to fall sleep at the beginning of treatment.
Supine Delay	30 sec	Time between detection of supine position and activation of stimulus.	Allows the patient to transition between left and right positions without applying any stimulus.
Stimulus Increase Period	60 sec	Time between consecutive increases of stimulus intensity.	Modifies the rate of stimulus increase for adapting it to each patient's needs.
Time at Maximum Power	20 min	Maximum time with the stimulus at 100% of intensity if patient does not change position.	Avoids excessive application of vibration stimulus.
Back to Sleep Delay	5 min	When the device detects the upright position and returns to any other position, it waits some minutes before applying again the stimulus.	Allows the patient to fall sleep again after getting up in the middle of the night.

3.3 VISUAL INDICATIONS AND WARNINGS

	TURN ON	
	RUNNING	
	TURN OFF	
	LOW BATTERY	Connect the device to the mains with the battery charger.
	CHARGING	The light turns OFF when the device is fully charged.
	UNSTABLE CLIPPING	Adhesive incorrectly clipped to the device. Check the clipping of both snap buttons to avoid device uncoupling.
	SEFT-TEST ERROR	Internal device error. Contact SIBEL's After Sales Service or your dealer.

 (vibration)

3.4 TRAVELLING WITH THE SORNIBEL

 Do not throw the packaging after using the device. It can be used for device transportation.

If you need to travel, you can bring your **Sornibel** with you for continuing your treatment. **Sornibel** can be easily and safely transported in your personal luggage by placing the device and accessories (battery charger, fastening adhesive box, etc.) into the original packaging.

4. CLEANING INSTRUCTIONS

The **Sornibel** device and the battery charger should be gently cleaned, with a cloth moistened with soapy (neutral) water or with 96° alcohol. Then they can be wiped dry. Particular care must be taken to ensure that no liquid enters the interior of the device or the connectors and connections. Abrasive substances or solvents must be avoided.

Fastening adhesives should not be cleaned and must be disposed after each use.

5. TECHNICAL SPECIFICATIONS

	Somnibel
Useful life	3 years
Size (L x W x H)	51.9 x 31.7 x 14.5 ± 0.5 mm
Weight	16.8 ± 1.0 gr
Power supply	Internal rechargeable Lithium-Polymer battery (> 3 nights @ 10h/night)
Enclosure protection degree:	IP22
Environmental conditions	Use: Temperature range: 5°C to 40°C (Note 1) / Relative humidity range: 35% to 85%(non-condensing) / Atmospheric pressure range: 700 hPa to 1060 hPa (from 2950 m to 350 m of altitude, approximately) Transport and storage conditions: Device (Note 2): Temperature range: -20°C to 60°C. / Humidity range: <93%(non-condensing). Adhesives: Temperature range: -25°C to +70°C / Humidity range: <93% (non-condensing).
Intensity Level at Maximum Stimulus	< 23 dB HL
A-weighted sound pressure level, L_{pA} A weighted sound power level, L_{WA}	DECLARED SINGLE-NUMBER NOISE EMISSION VALUES in accordance to ISO 4871 (Note 3) < 36 dBA < 35 dBA
USB 2.0 connectivity	Only for battery charge
Automatic ON/OFF	X
Initial Self-Test	X
Operational modes	Therapy
Monitored signals	Sleeping position
Compliance Enhancement Settings (Note 4)	

	Somnibel
Auto-adjusting Stimulus Levels	0%, 25%, 50%, 75%, 100% (< 110 dB)
EasySleep Delay	15 min
Supine Delay	30 secs
Stimulus Increase Period	60 secs
Time at Maximum Power	20 min.
Supine threshold angle	45°
Visual indications (Note 5)	Device Turn ON / Device Turn OFF / Device running / Low battery / Charging battery / Self-Test error Incorrect Clipping

Notes:

1. Maximum temperature is restricted to 37°C only when the device is used during the treatment (with the adhesive). Thus, the temperature on the adhesive will be below 40°C (maximum recommended temperature for the adhesive) and below 41°C (according to EN60601-1 standard). / **2.** Transport and storage conditions of the device are reduced with respect to the IEC 60601-1-11 standard due to the specifications of the internal lithium –polymer battery, vibration motor and switches. / **3.** Values measured at a distance of 0.25 m for the maximum intensity level with a test method equivalent to ISO 3744:2010. / **4.** See section COMPLIANCE ENHANCEMENT SETTINGS for definitions. / **5.** See section VISUAL INDICATIONS AND WARNINGS

6. APPLICABLE DIRECTIVES AND STANDARDS

- European directive concerning medical devices 93/42/EEC
- Rohs Directive: 2011/65/EU.
- Waste disposal according to WEEE Directive 2012/19/UE.
- Regulation EC 1272/2008 on classification, labeling and packaging of substances and mixtures (REACH).
- Packaging and packaging waste directive 94/62/ EC modified by 2004/12/EC
- Batteries waste directive 2006/66/EC modified by 2013/56/EU.

- Quality (EN ISO 13485:2016+AC:2018, EN ISO 9001:2015)
- Risk management (EN ISO 14971:2012)
- Safety Medical devices (EN 60601-1:2006+AC:2010+A1:2013+AC:2014)
- Safety Medical devices – Batteries (EN 62133:2013).
- Safety Medical devices – Home use (EN 60601-1-11:2015).
- Electro-magnetic Compatibility (EN 60601-1-2:2015; home and healthcare environment)
- Biocompatibility: Biological evaluation of medical devices (EN ISO 10993-1:2009+AC:2010)
- Usability (EN 60601-1-6:2010+A1:2015 and EN 62366:2008+A1:2015)
- Software of medical devices (EN 62304:2006+AC:2008+A1:2015)
- Documentation and information (EN 1041:2008 and EN ISO 15223-1:2016)

7. SYMBOLS

	MANUFACTURER (manufacture date, name and address of the manufacturer)		DIRECT CURRENT		ATMOSPHERIC PRESSURE LIMITATION
	DATE OF MANUFACTURE	IPN.N.	IP22 (SEE SECTION WARNINGS)		HUMIDITY LIMITATION
	CATALOGUE NUMBER		DEVICE UPRIGHT POSITION		KEEP DRY
	SERIAL NUMBER		BATCH CODE		FRAGILE, HANDLE WITH CARE
CE 0197	CE MARKING, NOTIFIED BODY NUMBER		USE BY DATE		RECYCLABLE CARDBOARD
	BF APPLIED PARTS		DOES NOT CONTAIN NATURAL RUBBER LATEX		DO NOT REUSE
	CAUTION, ADDITIONAL WARNINGS IN ACCOMPANYING DOCUMENTS		WASTE OF ELECTRICAL AND ELECTRONIC EQUIPMENT		THIS WAY UP PACKAGING
	WARNING, RISK IDENTIFIED		TEMPERATURE LIMITS		
	CONSULT INSTRUCTIONS FOR USE				

8. BIBLIOGRAPHY

1. Review by Veasey SC, Guilleminault C, Strohl KP et al. Medical Therapy for Obstructive Sleep Apnea: A Review by the Medical Therapy for Obstructive Sleep Apnea Task Force of the Standards of Practice Committee of the American Academy of Sleep Medicine. *Sleep* 2006;29(8):1036-1044.
2. Penzel T, Möller M, Becker HF, et al. Effect of sleep position and sleep stage on the collapsibility of the upper airways in patients with sleep apnea. *Sleep* 2001;24:90-95.
3. Oksenberg A, Khamaysi Y, Silverberg DS, et al. Association of body position with severity of apneic events in patients with severe nonpositional obstructive sleep apnea. *Chest* 2000;118:1018-1024.
4. Sahlin C, Svanborg E, Stenlund H, et al. Cheyne-Stokes respiration and supine dependency. *Eur Respir J* 2005; 25: 829-833.
5. Cartwright RD. Effect of sleep position on sleep apnea severity. *Sleep* 1984; 7: 110-14.
6. Oksenberg A, Silverberg DS, Arons E, et al. Positional vs nonpositional obstructive sleep apnea patients. Anthropomorphic, nocturnal polysomnographic and multiple sleep latency test data. *Chest* 1997;112:629-639.
7. Mador MJ, Kufel TJ, Magalang UJ, et al. Prevalence of positional sleep apnea in patients undergoing polysomnography. *Chest* 2005;128:2130-2137.
8. Nakano H, Ikeda T, Hayashi M, et al. Effects of Body Position on Snoring in Apneic and Nonapneic Snorers. *Sleep* 2003; 2:169-172.
9. Skinner MA, Kingshott RN, Filsell S, et al. Efficacy of the 'tennis ball technique' versus nCPAP in the management of position-dependent obstructive sleep apnoea syndrome. *Respirology* 2008; 13: 708-715.
10. Epstein LJ, Kristo D, Strollo PJ, et al. Clinical guideline for the evaluation, management and long-term care of obstructive sleep apnea in adults. *J Clin Sleep Med* 2009;5(3):263-276.
11. Hidalgo Armas L, Turino C, Cordero-Guevara J et al. A new postural device for the treatment of positional obstructive sleep apnea. A pilot study. *Respir Med* 2019; 151:111-117.

Annex 1: ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacturer's declaration – electromagnetic emissions		
Sornibel is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - Guidance
RF (Radiated) emissions CISPR 11 (EN 55011)	Group 1 Class B	Sornibel uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF (Conducted) emissions CISPR 11 (EN 55011)	Grupo 1 Clase B	Sornibel works with batteries. Applicable to the charger.
Harmonic emissions EN-IEC 61000-3-2	Not applicable	Sornibel works with batteries. Power of the charger is lower than 75W.
Voltage fluctuations / Flicker emissions EN-IEC 61000-3-3	Yes	Sornibel works with batteries. Applicable to the charger.

Guidance and manufacturer's declaration – electromagnetic emissions

Somnibel is intended for use in the electromagnetic environment specified below. The customer or the user of **Somnibel** should assure that it is used in such an environment. Tests EN-IEC 61000-4-4 y -4-5 are applicable to AC/DC power inputs, and input/output signal. Test EN-IEC 61000-4-6 additionally applies to patient connections. Test EN-IEC 61000-4-11 only applies to AC power inputs.

Immunity test	EN-IEC 60601 test level	Compliance level	Electromagnetic environment –Guidance
Electrical fast transient/burst EN-IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines Not applicable	Somnibel works with batteries. The input/output line cables are shorter than 3 meters long.
Surge EN-IEC 61000-4-5	±0.5, ±1 kV differential ±0.5, ±1, ±2 kV common mode	±0.5, ±1 kV differential ±0.5, ±1, ±2 kV common mode	Somnibel works with batteries.
Conducted RF EN-IEC 61000-4-6	3 Vrms de 150KHz a 80 MHz 6 Vrms in ISM bands 80% AM a 1kHz	3 Vrms 6 Vrms	
Voltage dips, short interruptions and voltage variations on power supply input lines EN-IEC 61000-4-11	<5 % Ut (>95 % dip in Ut) for 0.5 cycle 40 % Ut (60 % dip in Ut) for 5 cycles 70 % Ut (30 % dip in Ut) for 25 cycles <95 % Ut (>5 % dip in Ut) for 5 seconds	Not applicable	Somnibel works with batteries.

NOTE Ut is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

Somnibel is intended for use in the electromagnetic environment specified below. The customer or the user of Somnibel should assure that it is used in such an environment. Below mentioned tests are applicable to enclosure. Additionally, test EN-IEC 61000-4-2 is applicable to input/output signal and patient connections.

Immunity test	EN-IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN-IEC 61000-4-2	±8 kV contact ±2, ±4, ±8, ±15 kV, air	±8 kV contact ±2, ±4, ±8, ±15 kV, air	
Radiated RF EN-IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	10 V/m	
Proximity fields from RF communications EN-IEC 61000-4-3	(View following table).		Portable and mobile RF communications equipment should be used no closer than 30cm to any part of Somnibel including cables. The distance is calculated according to: $d = \left[\frac{6}{E} \right] \sqrt{P}$ where P is the output power rating of the transmitter in watts (W), d is the recommended separation in meters (m) and E is the immunity level in V/m.
Power frequency (50 / 60 Hz) magnetic field EN-IEC 61000-4-8	30 A/m	30 A/m	Only applicable to devices sensitive to magnetic fields.

Proximity fields EN-IEC 61000-4-3					
Frequency (MHz)	Band (MHz)	Modulation	Power (W)	Distance (m)	E (V/m)
385	380-390	Pulse 18Hz	1.8	0.3	27
450	430-470	FM D: +/-5kHz 1kHz sinus	2	0.3	28
710 745 780	704-787	Pulse 217Hz	0.2	0.3	9
810 870 930	800-960	Pulse 18Hz	2	0.3	28
1720 1845 1970	1700-1990	Pulse 217Hz	2	0.3	28
2450	2400-2570	Pulse 217Hz	2	0.3	28
5240 5500 5785	5100-5800	Pulse 217Hz	0.2	0.3	9

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All rights reserved.