
WIRKPRINZIP / METHOD OF ACTION / PRINCIPIO DE ACCIÓN / PRINCIPIO DI FUNZIONAMENTO / PRINCIPE DE FONCTIONNEMENT

RC-FIT® CLASSIC

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RC-FIT® CLASSIC

A) TRAINING & THERAPY WITH VIBRATIONS

The same principle applies to sports and everyday life: whoever has more oxygen available has a greater ability to function. The inlay with integrated valve produces oscillations that are acoustically perceived as sounds.

B) ANCHOR FOR SHORTNESS OF BREATH

As with the pursed lips technique, the RC-FIT® CLASSIC can stabilise the bronchi in cases of shortness of breath, through a deliberately slow exhalation. Hyperventilation can be reduced or prevented.


USE

The use of the RC-FIT® CLASSIC is not dependent on the force of gravity. The device may thus be used in a reclining, sitting, standing or moving position.

- I. Please inspect the device for foreign bodies or damages prior to every use.
II. Always start your first training session in the lowest intensity level "MINIMUM" (largest-possible triangle).

Table with columns: Average pressure during exhalation/Inhalation intensity level and four columns of flow rate values.

III. A) Place the RC-FIT® CLASSIC in your mouth and close your lips around the mouthpiece.

- A) If the sound at the start of each breathing manoeuvre is loud and then quickly becomes quieter, please adapt your breathing accordingly.
B) Make sure that the sound when breathing through the RC-FIT® CLASSIC is uniform and continuous.
C) Depending on your training condition, take a 1-minute break after 1½ minutes.

If the training is too light, or the sound is too loud even after 15 breaths, increase the intensity (resistance) by rotating the mouthpiece, see Fig. II.

CLEANING

- 1. Cleaning: For cleaning purposes, pull off the mouthpiece from the intensity regulator.
2. Drying: In order to prevent mould formation, all parts must be completely dried.
3. Assembly: Before assembly, please wash your hands and ensure that all parts of the RC-FIT® CLASSIC are clean.

FREQUENCY OF USE

Recommended therapy intervals for the respiration: RC-FIT® CLASSIC: 2 x daily 15 breaths, 7 days a week. Additional use when you notice heavy mucus obstruction or occasional shortness of breath.

- TIPPS:
Integrate the therapy and training into your daily routine, especially when you notice very little mucus, e.g. always before brushing your teeth.
To obtain a better effect from your medication spray, integrate this inhalation into your therapy with the RC-FIT® CLASSIC.

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